

Two-dimensional versus One-dimensional Transient Elastography: Benefits of Ultrasound Imaging Based Processing for Liver Stiffness Measurements

A. BESSON¹, B. HERIARD-DUBREUIL¹, V. DE LEDINGHEN², D. COHEN-DUTARTRE³, F. MANON², J. ABIVEN², A.-L. DE ARAUJO², R. HOUIMADI⁴, J. DUPUY⁴, J. FOUCHER², J. GAY¹, C. COHEN-BACRIE¹

1

Introduction

- Fibroscan® VCTE™: **1D** tracking of shear wave displacement generated with 50 Hz vibration + **A/M-Modes**
- Hepatoscope™ 2DTE: **2D** tracking of shear wave displacement generated with 50 Hz vibration + **B-Mode imaging**

2

Aim

To highlight benefits of 2DTE vs. 1D TE

- Shear wave propagation measured in the **direction of shear wave propagation** vs. in the **direction of ultrasound beam**
- **Quality control** of shear wave front
- **Visualization of the liver** to confirm the location of the shear wave measurement

3

Method

- Prospective single center study (NCT04782050)
- 96 patients referred to routine outpatient hepatology consultation for chronic liver disease
- 4 Hepatoscope™ liver exams
- 1 expert and 1 novice operators blinded to any median value
- 1D TE (LSM_{1D}) and 2DTE (LSM_{2D}) reprocessing
- Linear regression analysis LSM_{1D} vs. LSM_{2D}
- Repeatability / reproducibility (ICC)

4

Results

		LSM_{1D}	LSM_{2D}
Inter-operator Reproducibility	ICC [CI 95%]	0.59 [0.33-0.68]	0.74 [0.60-0.84]
Intra-operator Repeatability (experts)	ICC [CI 95%]	0.53 [0.33-0.68]	0.84 [0.73-0.91]




Figure 1: From left to right. Shear-Wave Propagation overlaid on B-Mode image as seen with 2DTE (leftmost image) and 1D TE (rightmost image). As the direction of propagation is not aligned with the analysis direction, 1D TE overestimates the shear-wave velocity.

5

Conclusions

2DTE has demonstrated its ability to **eliminate average over-estimation bias** of LSM by 1D TE

2DTE has led to **increased intra- and inter-operator reproducibility**

Further study against histology needed to demonstrate that 2DTE **reduces false positives of 1D TE**

6

Contact information

Adrien Besson, adrien.besson@e-scopics.com

2DTE: Seeing Is Believing