




GOOD DIAGNOSTIC PERFORMANCE OF HEPATOSCOPE™ IN DIABETOLOGY CLINICS FOR THE SCREENING OF MASLD PATIENTS WITH HIGH RISK OF ADVANCED LIVER FIBROSIS

C. Caussy^{1,2,3}, B. Vergès⁴, F. Subtil⁵, A. Abichou-Klich⁵, V. Hervieu⁶, L. Milot⁷, B. Ségrestin¹, S. Bin⁸, A. Rouland⁴, D. Delaunay¹, S. Hadjadj⁹, B. Cariou⁹, C. Primot⁹, J-M Petit⁴, P. Moulin^{2,10}, S. Charrière^{2,10}, E. Disse^{1,2}

1.Hospices Civils de Lyon, Département Endocrinologie, Diabète et Nutrition, Hôpital Lyon Sud, Pierre-Bénite, France; 2. Univ Lyon, CarMen Laboratory, INSERM, INRA, INSA Lyon, Université Claude Bernard Lyon 1, Pierre-Bénite, France ; 3. IHU Lyon - Institut d'Hépatologie de Lyon, France.; 4. Department of Endocrinology, Diabetes and Metabolic Disorders, Dijon University Hospital, INSERM Unit, LNC-UMR 1231, University of Burgundy, Dijon, France. 5.Hospices Civils de Lyon,Service de Biostatistique, Lyon, France; Université Claude Bernard Lyon 1, UMR 5558, CNRS, Villeurbanne, France 6. Hospices Civils de Lyon, Biopathology of Tumours, GH Est Hospital, Bron, France. 7. Hospices Civils de Lyon, Service de Radiologie, Hôpital Edouard Herriot, 5 Place d'Arsonval, Lyon, France; 8. Hospices Civils de Lyon, Pôle de Santé Publique, Service Recherche et Épidémiologie Cliniques, Lyon, France 9. Nantes Université, CHU Nantes, CNRS, INSERM, l'institut du thorax, Nantes, France; 10. Hospices Civils de Lyon, Fédération d'Endocrinologie, Diabète et Nutrition, Hôpital Cardiovasculaire, 69500 Bron, France



Contact information
cyrielle.caussy@chu-lyon.fr

Introduction

The screening for metabolic dysfunction-associated steatotic liver diseases (MASLD)-related advanced fibrosis (AF) is recommended in patients with type 2 diabetes (T2D) and obesity with metabolic co-morbidities. The Hepatoscope® is an ultraportable ultrasound system with two-dimensional transient elastography (2DTE) for liver stiffness measurement (LSM) (1-2).

Aim


To evaluate the diagnostic performance of a new ultrasound imaging point-of-care device: Hepatoscope for the detection of patients with high risk of AF in diabetology.

Method

- Ancillary study from a prospective multicenter NAFLD-Care study (NCT04435054) (3) including participants willing to participate in the sub-study for the assessment of the diagnostic performance of the Hepatoscope™
 - Patients with MASLD and T2M and/or obesity, age 40-80 years, BMI < 40 kg/m².
 - 4 diabetes clinics in France
 - Enrollment between December 2022 and September 2024
 - Exclusion of all other liver etiologies
- LSM was performed on the same day using 2 ultrasound devices at the fasting state using :
 - Vibration controlled transient elastography (VCTE™) by FibroScan® and
 - Hepatoscope: Prototype version 2DTE-V1 for acquisition / post-acquisition reprocessing of data with CE-certified version (2DTE-V2)
- High risk of AF was determined using a hierarchical composite criteria: 1. Liver biopsy ≥ F3 if available or 2. Magnetic resonance elastography (MRE) ≥ 3.62 kPa or 3. VCTE ≥ 12 kPa as previously published (3)
- Statistical analysis: Spearman's correlation between VCTE and 2DTE-V1 or 2DTE-V2 in participants with valid VCTE (IQR/M ≤ 0.30) and AUROC and 95%CI for the detection of AF.
- Exploratory subgroup analysis performed in 48 participants with magnetic resonance imaging proton density fat fraction (MRI-PDFF) and backscattering coefficient (BSC) by Hepatoscope™ for the quantification of hepatic fat content.

Results

1.Study flowchart



2.Baseline characteristics of the study populations

	Population V1 (n=296)	Population V2 (n=119)	p-value V1 vs V2
Age (years)	61.0 [53.2 - 67.0]	60.0 [53.0 - 66.0]	0.40
Female (n, %)	127 (42.9)	44 (37.0)	0.32
T2D (n, %)	247 (83.4)	96 (80.7)	0.57
BMI (kg/m²)	32.5 [29.2 - 67.0]	32.2 [29.0 - 35.8]	0.54
Obesity (n, %)	210 (70.9)	82 (69.7)	0.81
Non-invasive tests			
FIB-4	1.34 [0.93 - 1.72]	1.32 [0.99 - 1.57]	0.64
LSM _{VCTE} (kPa)	6.1 [4.8 - 8.0]	6.8 [5.0 - 8.6]	0.14
FibroScan XL Probe (n,%)	126 (42.6)	36 (30.0)	0.0011
LSM _{MRE} (kPa) n=59	2.64 [2.28 - 3.00]	2.69 [2.31 - 3.36]	0.18
LSM _{2DTE-V1} (kPa)	5.7 [4.4 - 7.5]	NA	NA
LSM _{2DTE-V2} (kPa)	NA	7.1 [4.4 - 7.5]	NA
MRI-PDFF (%) n=59	10.0 [6.5 - 23.0]	10.0 [6.9 - 23.2]	0.61
CAP (dB/m)	306 [273 - 343]	303 [266 - 335]	0.41
BSC (dB/m.sr)	NA	-19.8 [-28.3 - -23.2]	NA
LSM_{VCTE} ≥ 8 (kPa) (n,%)	74 (25.0)	38 (31.7)	0.18
Presence of AF (n,%)*	70 (8.9)**	15 (12.7)	0.23

* Presence of AF determined using hierachical composite criteria as follow: **n=25 with liver biopsy with n=10 ≥ F3; n=51 MRI with n=7 ≥ 3.6 kPa or overt cirrhosis; n=220 VCTE (n=14≥ 12kPa).

3. Comparaison between LSM_{2DTE} and LSM_{VCTE}




Figure 1. Comparison of LSM_{2DTE} and LSM_{VCTE}.
A. Significant Spearman's correlation between LSM_{2DTE-V1} and LSM_{VCTE}
B. Improved correlation between LSM_{2DTE-V2} and LSM_{VCTE}
C. No significant bias between LSM_{2DTE-V2} and LSM_{VCTE} as shown by the Bland-Altman analysis.

5. Hepatoscope performance to risk stratify patient for advanced fibrosis




Figure 2. AUROC of Hepatoscope 2DTE-V2
A. to identify patients with LSM_{VCTE} > 8 kPa
B. to identify patients with high risk for advanced fibrosis (composite endpoint)

Conclusions

- This proof-of-concept study indicates that Hepatoscope™ 2DTE has a good diagnostic performance for the identification of patients with high risk of AF among patients with T2D and/or obesity screened in diabetology.
- Further studies in larger cohort using the 2DTE-V2 by Hepatoscope™ are needed to determine the optimal cut-offs.

References

- Hériard-Dubreuil B et al. Ultraportable Quantitative Ultrasound for Hepatic Steatosis Assessment. Ultrasound Med Biol. 2024 Dec;50(12):1842-1848. doi: 10.1016/j.ultrasmedbio.2024.08.008.
- Serdjebi C. Hepatoscope 2DTE's image-based quality index enhances applicability and repeatability of liver stiffness measurement. Clin Res Hepatol Gastroenterol. 2025 Feb;49(2):102532. doi: 10.1016/j.clinre.2025.102532.
- Caussy C et al. Screening for Metabolic Dysfunction-Associated Steatotic Liver Disease-Related Advanced Fibrosis in Diabetology: A Prospective Multicenter Study. Diabetes Care. 2025 Jan 31

4. Risk stratification for advanced fibrosis





Figure 2. Composite score for stratification for risk of advanced fibrosis, based on liver biopsy, MRE and LSM_{VCTE}.

6. Risk stratification of patient for advanced fibrosis using LSM

	LSM < 8 kPa	8 ≤ LSM < 12	12 kPa ≤ LSM	p-value
LSM _{VCTE}	81	27	11	ns
LSM _{2DTE}	73	31	15	

Table 1. Risk stratification for advanced fibrosis based on LSM.
No significant change of distribution (Fisher's exact test).

7.Correlation between BSC and MRI-PDFF



Acknowledgements

The study was funded with the support of E-Scopics.

